



Malyugin Ring/ I-Ring: Scrolls/ Channels/ Pockets

Brilliant !!! ... Yet the biggest Hurdle !!



Radically Different Concept - The Inventive Step!


Malyugin Ring/ I - Ring



Scrolls/ Channels/ Pockets

- **Biplanar Structure**
 - Snags Incision, Requires Injector, Thick Profile
- Scroll straddles Pupil margin

B-HEX (Bhattacharjee Ring)



Notches & Flanges

- **Single Plane**
 - Injector Obviated, Thin Profile
- Pupil margin straddles the Notch



Iris can be bent without harm